

Anti-BLyS/TNFSF13B antibody



Catalog Number: 106634

Product name

Anti-BLyS/TNFSF13B antibody

Immunogen

[Human BLyS/TNFSF13B recombinant protein](#)

Specificity

Human BLyS / TNFSF13B

Antibody description

Rabbit monoclonal to BLyS/TNFSF13B

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human BLyS / TNFSF13B (rh BLyS / TNFSF13B; Q9Y275-1; Ala134-Leu285) and conjugated with PE under optimum conditions, the unreacted PE was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM

Validations



BLyS / TNFSF13B Antibody (PE), Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Human BAFF(CD257) expression on U937 cells. Cells were stained with PE-conjugated anti-Human BAFF(CD257). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.