

Anti-TBC1D3C antibody



Catalog Number: 176334

Product name

Anti-TBC1D3C antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TBC1D3C

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tbc1d3. aa range:481-530

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using TBC1D3A/B/C Polyclonal Antibody