

Anti-TBC1D3C antibody



Catalog Number: 176334

Product name

Polyclonal

Anti-TBC1D3C antibody

Ig Type

Specificity

Rabbit IgG

Human

Applications

Antibody description

WB, IF

Rabbit polyclonal antibody to TBC1D3C

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tbc1d3. aa range:481-530

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1.; Western Blot analysis of various cells using TBC1D3A/B/C Polyclonal Antibody