

# NOB1 antibody

Catalog Number: 107334

## Product name

NOB1 antibody

## Specificity

Human; other species not tested.

## Antibody description

NOB1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, human placenta tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 47 kDa

## Preparation

This antibody was obtained by immunization of NOB1 recombinant protein (Accession Number: BC000050). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107334(NOB1 Antibody) at dilution of 1:3000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107334(NOB1 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107334(NOB1 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9).

# NOB1 antibody



Catalog Number: 107334

---

No:107334(NOB1 Antibody) at dilution of 1:200

(under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).