

Human XPNPEP1/XPNPEPL/XPNPEPL1 (His tag) recombinant protein



Catalog Number: 518334

General Information

Gene Name Synonym

Xaa-Pro Aminopeptidase 1, Aminoacylproline Aminopeptidase, Cytosolic Aminopeptidase P, Soluble Aminopeptidase P, sAmp, X-Pro Aminopeptidase 1, X-Prolyl Aminopeptidase 1 Soluble, XPNPEP1, XPNPEPL, XPNPEPL1

Protein Construction

A DNA sequence encoding the Human XPNPEP1/XPNPEPL/XPNPEPL1 Pro2-His623 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human XPNPEP1/XPNPEPL/XPNPEPL1 (His tag)
recombinant protein