

# IL2RG antibody



Catalog Number: 111734

## Product name

IL2RG antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

IL2RG Rabbit Polyclonal antibody. Positive FC detected in Jurkat cells. Positive WB detected in Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 42kd

## Preparation

This antibody was obtained by immunization of IL2RG recombinant protein (Accession Number: NM\_000206). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111734(IL2RG antibody) at dilution of 1:400



1X10<sup>6</sup> Jurkat cells were stained with 0.2ug IL2RG antibody (Catalog No:111734, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.