

# STRADB antibody



Catalog Number: 115734

---

## Product name

STRADB antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

STRADB Rabbit Polyclonal antibody. Positive WB detected in Raji cells, mouse heart tissue. Observed molecular weight by Western-blot: 28-31 kDa; 38-42 kDa

## Preparation

This antibody was obtained by immunization of STRADB recombinant protein (Accession Number: NM\_018571). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115734(STRADB Antibody) at dilution of 1:300