

# Anti-p40 Rabbit antibody



Catalog Number: 166734

## Product name

Anti-p40 Rabbit antibody

## Specificity

Human

## Antibody description

p40 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human p40

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/20

## Validations



Western blot detection of p40 in K562, HeLa cell lysates using p40 Rabbit pAb(1:1000 diluted). Predicted band size: 41 kDa. Observed band size: 40 kDa.