

VDR antibody



Catalog Number: 116734

Product name

Polyclonal

VDR antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

VDR Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 55 kDa

Dilutions


Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of VDR recombinant protein (Accession Number: XM_047429500). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116734(VDR Antibody) at dilution of 1:1000

Clonality