

ACDP1 antibody



Catalog Number: 107734

Product name

ACDP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACDP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 100-110 kDa

Preparation

This antibody was obtained by immunization of ACDP1 recombinant protein (Accession Number: BC098103). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107734(CNNM1 Antibody) at dilution of 1:1000