

KIAA0141 antibody



Catalog Number: 112044

Product name

KIAA0141 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

KIAA0141 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of KIAA0141 recombinant protein (Accession Number: XM_005268548). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112044(KIAA0141 Antibody) at dilution of 1:300