

H2N2 Nucleoprotein (His Tag)



Catalog Number: 503044

General Information

Gene Name Synonym

Nucleocapsid protein

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Ann Arbor/6/1960 (H2N2)) nucleoprotein (AAM75159.1) (Met 1-Asn 498) was fused with a polyhistidine tag at the C-terminus.

Organism

H2N2

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met 1

Molecular Mass

The recombinant Influenza A virus (A/Ann

Arbor/6/1960(H2N2)) nucleoprotein comprises 508 amino acids and has a predicted molecular mass of 57 kDa. It migrates as an approximately 52 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 8.0, 20% gly, 1mM DTT

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H2N2 Nucleoprotein / NP Protein (His Tag) SDS-PAGE