

SAFB2 antibody



Catalog Number: 115044

Product name

SAFB2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SAFB2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, BxPC-3 cells, HeLa cells, MCF7 cells, rat pancreas tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of SAFB2 recombinant protein (Accession Number: BC025279). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115044(SAFB2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115044(SAFB2 antibody) at dilution of 1:50 (under 40x lens)