

Anti-PPP1R14C antibody



Catalog Number: 176044

Product name

Anti-PPP1R14C antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to PPP1R14C

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ppp1r14c. aa range:51-100

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using KEPI Polyclonal Antibody diluted at 1: 1000