

Anti-Caspase-9 Rabbit antibody



Catalog Number: 168044

Product name

Anti-Caspase-9 Rabbit antibody

Specificity

Human

Antibody description

Caspase-9 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Caspase-9

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC

Validations



Western blot detection of Caspase-9 in HeLa,A549,HL-60,U2OS using Caspase-9 Rabbit mAb(1:1000 diluted)