

# APRT antibody

Catalog Number: 108044

## Product name

APRT antibody

## Specificity

Human; other species not tested.

## Antibody description

APRT Rabbit Polyclonal antibody. Positive WB detected in A549 cells, A431 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of APRT recombinant protein (Accession Number: BC106894). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:108044(APRT antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed A549 cells using Catalog No:108044(APRT Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)