

Anti-S100A12 antibody

Catalog Number: 100444

Product name

Anti-S100A12 antibody

Immunogen

[Human S100A12 recombinant protein](#)

Specificity

Human S100A12

Antibody description

Rabbit monoclonal to S100A12

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human S100A12 (rh S100A12; NP_005612.1; Met1-Glu92).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 2-10 µg/mL

Validations



S100A12 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human S100A12 in human tonsil with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).