

# FABP7-Specific antibody

Catalog Number: 110444

## Product name

FABP7-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

FABP7-Specific Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in human brain tissue, U-251 cells. Observed molecular weight by Western-blot: 14-15 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001446). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110444(FABP7-Specific antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110444(FABP7-Specific antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110444(FABP7-Specific antibody) at dilution of

# FABP7-Specific antibody



Catalog Number: 110444

---

1:100 (under 40x lens)