

# HP1BP3 antibody



Catalog Number: 111444

## Product name

HP1BP3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HP1BP3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of HP1BP3 recombinant protein (Accession Number: BC045660). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111444(HP1BP3 Antibody) at dilution of 1:600