

Phospho-PRKCZ-T410 pAb



Catalog Number: 166444

Product name

Phospho-PRKCZ-T410 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-PRKCZ-T410

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding T410 of human PRKCZ

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions

IHC 1:50 - 1:100

Validations



Immunohistochemistry - Phospho-PRKCZ-T410 pAb

Immunohistochemistry of paraffin-embedded human lung carcinoma tissue, using Phospho-PRKCZ-T410 antibody .