

# Anti-DDIAS antibody

Catalog Number: 176444

## Product name

Anti-DDIAS antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to DDIAS

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human noxin. aa range:111-160

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.; Western Blot analysis of 823, 293T, AD293, Hela, HepG2 cells using NOXIN Polyclonal Antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000