

Anti-GRK3 Rabbit antibody



Catalog Number: 167444

Product name

Anti-GRK3 Rabbit antibody

Specificity

Human

Antibody description

GRK3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human GRK3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000-1/5000

IHC: 1/20-1/100

ICC/IF: 1/20-1/50

Validations



Western blot detection of GRK3 in Jurkat cell lysates using GRK3 Rabbit pAb(1:1000 diluted).Predicted band size:80KDa.Observed band size:80KDa.