

BEND4 antibody

Catalog Number: 108444

Product name

BEND4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BEND4 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human heart tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 55-58 kDa; 33-36 kDa

Preparation

This antibody was obtained by immunization of BEND4 recombinant protein (Accession Number: BC128592). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:117117(BEND4 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:117117(BEND4 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:117117(BEND4 Antibody) at dilution of 1:50

BEND4 antibody



Catalog Number: 108444

(under 40x lens)