

# COMMD5 antibody

Catalog Number: 109444

## Product name

COMMD5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

COMMD5 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse kidney tissue, HEK-293 cells, human heart tissue, human stomach tissue, mouse heart tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 25-28 kDa

## Preparation

This antibody was obtained by immunization of COMMD5 recombinant protein (Accession Number: NM\_001081004). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109444(COMMD5 antibody) at dilution of



IP Result of anti-COMMD5 (IP:Catalog No:109444, 4ug; Detection:Catalog No:109444 1:500) with mouse heart tissue lysate 3200ug.