

# Decorin antibody



Catalog Number: 109844

## Product name

Decorin antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Decorin Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, MCF7 cells. Positive IHC detected in human ovary tissue. Observed molecular weight by Western-blot: 45-48 kDa

## Preparation

This antibody was obtained by immunization of Decorin recombinant protein (Accession Number: NM\_001920). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109844(DCN Antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human ovary using Catalog No:109844(DCN antibody) at dilution of 1:50 (under 40x lens)