

DYRK1B antibody



Catalog Number: 110144

Product name

DYRK1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DYRK1B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HeLa cells, human colon tissue, human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of DYRK1B recombinant protein (Accession Number: XM_005259398). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110144(DYRK1B antibody) at dilution of 1:500