

# Phospho-MYL9-S19 pAb



Catalog Number: 166544

## Product name

Phospho-MYL9-S19 pAb

## Specificity

Rat

## Antibody description

Polyclonal antibody to Phospho-MYL9-S19

## Preparation

Antigen: A synthetic phosphorylated peptide around S19 of human MYL9 (NP\_006088.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - Phospho-MYL9-S19 pAb

Western blot analysis of extracts of C6 cells, using Phospho-MYL9-S19 antibody at 1:1000 dilution. C6 cells were treated by CIP(20uL/400ul) at 37u2103 for 1 hour. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 90s.