

# WTAP antibody

Catalog Number: 107544

## Product name

WTAP antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

## Specificity


Human, Mouse; other species not tested.

IF: 1:20-1:200

## Validations

### Antibody description

WTAP Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells, human testis tissue, Jurkat cells, mouse testis tissue. Positive IF detected in A431 cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 50 kDa




### Preparation

This antibody was obtained by immunization of WTAP recombinant protein (Accession Number: NM\_152858). Purification method: Protein A purified.

Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107544(WTAP antibody) at dilution of 1:50 (under 10x lens)

### Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



### Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:107544(WTAP antibody) at dilution of 1:50

### Clonality

Monoclonal

### Ig Type

Mouse IgG1


### Applications

ELISA, IHC, IF, WB

### Dilutions

Recommended Dilution:

(under 40x lens)




# WTAP antibody

Catalog Number: 107544

---

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107544(WTAP Antibody) at dilution of 1:2000



Immunofluorescent analysis of A431 cells using Catalog No:107544(WTAP Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)