

Anti-IL12A / IL-12A antibody



Catalog Number: 100944

Product name

Anti-IL12A / IL-12A antibody

Specificity

Mouse IL12A

Antibody description

Rabbit polyclonal to IL12A / IL-12A

Preparation

Produced in rabbits immunized with purified, a synthetic peptide corresponding to the N-terminus of the mouse IL12A . Mouse IL12A specific IgG was purified by Mouse IL12A affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-10 µg/ml

Validations



Lanes	A	B	C	D
Items				
Sample (whole cell lysate)	NIH-3T3	Hu-T78	Raw264.7	MDA-MB-231
Sample Volume (µg/lane)	30	30	30	30
Gel	13%SDS-PAGE reducing gel			
Recommended Concentration	5-10µg/ml			
Secondary Antibody	Dylight 800-labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution			
	Developed using Odyssey imaging system.			
Explanation	Predicted band size : 25 kDa Observed band size : 28 kDa Additional bands at : 40 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).			

IL12A Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot