

# PLCB3 antibody



Catalog Number: 113944

## Product name

PLCB3 antibody

## Specificity

Human; other species not tested.

## Antibody description

PLCB3 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A431 cells, COLO 320 cells, HEK-293T cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 150 kDa

## Preparation

This antibody was obtained by immunization of PLCB3 recombinant protein (Accession Number: NM\_000932). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113944(PLCB3 Antibody) at dilution of 1:600