

Anti-TPD52L2 antibody

Catalog Number: 175944

Product name

Anti-TPD52L2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TPD52L2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tpd54. aa range:141-190

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using D54 Polyclonal Antibody