

Anti-Kappa Opioid Receptor Rabbit antibody



Catalog Number: 168944

Product name

Anti-Kappa Opioid Receptor Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Rat

Applications

WB, IHC-P, IP

Antibody description

Kappa Opioid Receptor Rabbit monoclonal antibody

Dilutions

1:1000

1:20

1:20

Preparation

Antigen: A synthetic peptide of mouse Kappa Opioid Receptor

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of Kappa Opioid Receptor in Rat Brain lysates using Kappa Opioid Receptor Rabbit mAb(1:1000 diluted). Predicted band size:43kDa.Observed band size:60kDa.