

# Anti-STRA8 antibody



Catalog Number: 142244

## Product name

Anti-STRA8 antibody

## Specificity

Human STRA8

## Antibody description

Rabbit Polyclonal to Human STRA8

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human STRA8, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 5-10 µg/ml

## Validations



Items	Lanes	A	B
Sample (whole cell lysate)		HEK293	293
Sample Volume (µg/lane)		30	30
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	5-10 µg/ml		
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.		
Explanation	Developed using Odyssey imaging system. Predicted band size : 37 kDa Observed band size : 35 kDa Additional bands at : 60 kDa and 40 kDa (We are unsure as to the identity of these extra bands.)		

Human STRA8 Western blot (WB) 14560