

PTHLH Polyclonal Antibody



Catalog Number: 163244

Product name

PTHLH Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PTHLH

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-175 of human PTHLH (NP_945315.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - PTHLH Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using PTHLH antibody at dilution of 1:200



(40x lens).

Immunohistochemistry - PTHLH Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using PTHLH antibody at dilution of



1:200 (40x lens).

Immunohistochemistry - PTHLH Polyclonal Antibody

PTHLH Polyclonal Antibody



Catalog Number: 163244

Immunohistochemistry of paraffin-embedded

human gastric cancer using PTHLH antibody at dilution of 1:200 (40x lens).