

# Human TLK2/PKU-ALPHA (aa 397-772) recombinant protein



Catalog Number: 504244

## General Information

### Gene Name Synonym

HsHPK; PKU-alpha; Tousled-like kinase 2

### Protein Construction

A DNA sequence encoding the human TLK2 (Q86UE8-1) (Leu397-Asn772) was fused with two additional amino acids (Gly & Pro) at the N-terminus.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Activity

Kinase activity untested

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Gly

### Molecular Mass

The recombinant human TLK2 consists of 378 amino acids and has a calculated molecular mass of 43.6 kDa. The recombinant protein migrates as an approximately 44 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Supplied as sterile 20mM Tris, 500mM NaCl, 3mM DTT, 10% glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human TLK2 / PKU-ALPHA Protein (aa 397-772)  
SDS-PAGE