

# SIRT3 antibody

Catalog Number: 115244

---

**Product name**

SIRT3 antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

SIRT3 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human spleen tissue, mouse kidney tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 28 kDa

**Preparation**

This antibody was obtained by immunization of SIRT3 recombinant protein (Accession Number: NM\_012239). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:1000

**Validations**

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115244(SIRT3 Antibody) at dilution of 1:300