

Anti-VEGFR2/Flk-1/CD309/KDR antibody



Catalog Number: 106244

Product name

Anti-VEGFR2/Flk-1/CD309/KDR antibody

Immunogen

[Mouse VEGFR2/Flk-1/CD309/KDR \(His Tag\) recombinant protein](#)

Specificity

Mouse VEGFR2 / Flk-1 / CD309 / KDR

Antibody description

Rabbit monoclonal to VEGFR2/Flk-1/CD309/KDR

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse VEGFR2 / Flk-1 / CD309 / KDR (rM VEGFR2 / Flk-1 / CD309 / KDR; P35918-1; Met1-Glu762).

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

FCM

Dilutions

FCM: 0.5-2 μ g/Test

Validations



VEGFR2 / Flk-1 / CD309 / KDR Antibody, Rabbit MAb, Flow Cytometry

Flow cytometric analysis of Mouse KDR(CD309) expression on bEnd.3 cells. Cells were stained with purified anti-Mouse KDR(CD309), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.