

Anti-GLB1/Beta-galactosidase Rabbit antibody



Catalog Number: 167244

Product name

Anti-GLB1/Beta-galactosidase Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, FC

Antibody description

GLB1/Beta-galactosidase Rabbit polyclonal antibody

Dilutions

WB: 1/2000-1/10000

FC: 1/20

Preparation

Antigen: A synthetic peptide of human GLB1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Validations



Western blot detection of GLB1 in K562,C6,3T3,Hela cell lysates using GLB1 Rabbit pAb(1:1000 diluted).Predicted band size:76kDa.Observed band size:80kDa.