

Anti-FYB (2H9) Mouse antibody



Catalog Number: 168244

Product name

Anti-FYB (2H9) Mouse antibody

Specificity

Human

Antibody description

FYB (2H9) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human FYB protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/500

Validations



Western blot detection of FYB in 10,15 and 20ug Jurkat cell lysate using FYB mouse mAb (1:500 diluted). Predicted band size:120KDa. Observed band size:120KDa.