

ARF3 antibody

Catalog Number: 108244

Product name

ARF3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARF3 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 18-20 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001659). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108244(ARF3 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:108244(ARF3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:108244(ARF3 Antibody) at dilution of 1:50 (under 40x lens)