

CEACAM3-Specific antibody



Catalog Number: 109244

Product name

CEACAM3-Specific antibody

Specificity

Human; other species not tested.

Antibody description

CEACAM3-Specific Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human stomach cancer tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001815). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109244(CEACAM3-Specific antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109244(CEACAM3-Specific antibody) at dilution of 1:50 (under 40x lens)