

Anti-GOLGA5 antibody



Catalog Number: 176644

Product name

Anti-GOLGA5 antibody

Dilutions

WB: 1:1,000

Specificity

IHC-P: 1:100

Human

Antibody description

Mouse monoclonal antibody to GOLGA5

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

WB, IHC-P

Validations




Fig1: Western blot analysis on HeLa cell lysates using anti-GOLGA5 Mouse mAb (Cat. # 176644#).




Fig2: Immunohistochemical analysis of paraffin-embedded gastric cancer tissue using anti-GOLGA5 Mouse mAb (Cat. # 176644#).