

# URM1 antibody



Catalog Number: 116644

## Product name

URM1 antibody

## Immunogen

[Human URM1 Recombinant protein \(GST tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

URM1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A549 cells, HepG2 cells, mouse liver tissue. Positive IHC detected in human hepat°Cirrhosis. Observed molecular weight by Western-blot: 11kd

## Preparation

This antibody was obtained by immunization of URM1 recombinant protein (Accession Number: NM\_030914). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116644(URM1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116644(URM1 antibody) at dilution of 1:100 (under 10x lens)



# URM1 antibody



Catalog Number: 116644

---

Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog

No:116644(URM1 antibody) at dilution of 1:100 (under 40x lens)