

URM1 antibody



Catalog Number: 116644

Product name

URM1 antibody

Immunogen

[Human URM1 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

URM1 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, A549 cells, HepG2 cells, mouse liver tissue. Positive IHC detected in human hepatocirrhosis. Observed molecular weight by Western-blot: 11kd

Preparation

This antibody was obtained by immunization of URM1 recombinant protein (Accession Number: NM_030914). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:116644(URM1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:116644(URM1 antibody) at dilution of 1:100 (under 10x lens)



URM1 antibody



Catalog Number: 116644

Immunohistochemical of paraffin-embedded
human hepatocirrhosis using Catalog

No:116644(URM1 antibody) at dilution of 1:100
(under 40x lens)