

# CXCR4 antibody



Catalog Number: 109644

## Product name

CXCR4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CXCR4 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in rat spleen tissue, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 60-70kd

## Preparation

This antibody was obtained by immunization of CXCR4 recombinant protein (Accession Number: NM\_003467). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



IP Result of anti-CXCR4 (IP:Catalog No:109644, 3ug; Detection:Catalog No:109644 1:500) with



HeLa cells lysate 1000ug.

rat spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109644(CXCR4 Antibody) at dilution of 1:1000