

Oligophrenin 1 antibody



Catalog Number: 113344

Product name

Oligophrenin 1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Oligophrenin 1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 92 kDa

Preparation

This antibody was obtained by immunization of Oligophrenin 1 recombinant protein (Accession Number: BC059393). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113344(OPHN1 antibody) at dilution of 1:800



IP Result of anti-OPHN1 (IP:Catalog No:113344, 3ug; Detection:Catalog No:113344 1:200) with mouse brain tissue lysate 8000ug.