

SLC5A11 antibody



Catalog Number: 115344

Product name

SLC5A11 antibody

Specificity

Human; other species not tested.

Antibody description

SLC5A11 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of SLC5A11 recombinant protein (Accession Number: NM_001352259). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115344(SLC5A11 Antibody) at dilution of 1:300