

# Anti-GTF3C3 antibody

Catalog Number: 176344

## Product name

Anti-GTF3C3 antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to GTF3C3

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tf3c3. aa range:101-150

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1.; Western Blot analysis of various cells using TFIIC102 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).