

# Anti-NKX2.2 (1C3) Mouse antibody



Catalog Number: 168344

## Product name

Anti-NKX2.2 (1C3) Mouse antibody

WB: 1/500 - 1/2000

FC: 1/200 - 1/400

## Specificity

Human

ELISA: 1/10000

## Antibody description

NKX2.2 (1C3) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human NKX2.2 expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type


Mouse IgG1

## Applications


WB, FC

## Dilutions

## Validations



Western blot analysis of extracts from HCT116 and MCF7 cell lysates using NKX2.2 mouse mAb (1:1000 diluted). Predicted band size: 30 kDa. Observed band size: 38 kDa.



Flow cytometric analysis of MCF-7 cells using NKX2.2 mouse mAb (green) and negative control (red).