

# FOSL2 antibody



Catalog Number: 110744

---

## Product name

FOSL2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOSL2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 31 kDa

## Preparation

This antibody was obtained by immunization of FOSL2 recombinant protein (Accession Number: BC008899). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110744(FOSL2 Antibody) at dilution of 1:300